
### Trouble Shooting Scheren

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
## Trouble Shooting Scissors

Maintenance and re-processing of scissors

## Hinge-breaks! Rust attack in the lock and on the surface!

Quite often, scissors are sent back to us for unjustified COMPLAINT. Such scissors present the following faults:

- Rust in the instrument's lock, breaks in the instrument's lock and especially on the tips.

In most cases (> 95%) we can exclude faults of material or workmanship. In nearly all cases the defects occur due to faults in the hygiene chain. Below you will find our statement for those faults.

### Instrument breaks - tips

Generally, it is important for all instruments to put them softly down and not to throw them! All hinge instruments have to be cleaned and sterilized in open position!

Do not close them! This is very important for scissors. If they are closed, the instruments' tips are under contact pressure. In the scissors' construction, while hardening, the scissor blades are bent against each other and twisted around. So when assembly is finished, only the blades have contact at the intersection point. Due to the warming during sterilization, the scissors steel stretches. During this process, the contact pressure on the tips of the scissors is rising. If the scissors are open, this will not happen! Especially surgical scissors have very fine tips, and stress cracks can occur in closed position. Such cracks can even be in the micro range, so that the instruments will only break under strain!

When cleaning, it is always important for all instruments to put them softly down and not to throw them! Falling down, cutting something for which the scissors are inappropriate, and mechanical stress when the scissors are in contact with other instruments can lead to hairline cracks.

Different concentrations of chlorides are always found in drinking water. They cannot be eliminated by a water softening process. When water dries, chlorides will concentrate and cause pitting corrosion on the surface. In order to avoid this, the final rinsing step should be carried out with completely desalinated (demineralized) water. The pitting corrosion primarily leads to breaks on the tips!


**SCHWERT**  
SWORD · ESPADA  
EPEE · SPADA

## Trouble Shooting Scissors

Maintenance and re-processing of scissors

### Rust in the lock

If you ask what "stainless steel" and "corrosion-resistant high-grade steel" mean, it is often assumed that so-called stainless steel materials are extremely immune materials. The corrosion resistance of "STAINLESS STEELS" is based on their alloy structure which forms passivated surface layers that in turn protect steel. However, external influences can destroy these protective layers and cause severe damage to the instruments.

The scissors have to be cooled down and oiled prior to the function check to avoid metallic abrasion. Metallic abrasion can cause a jamming of hinges and lead to corrosion.

This is often a creeping process, primarily causing breakage of the scissors in the lock!

Prior to the function check, the hinges of scissors (or forceps or needle holders) have to be oiled with a care product based on paraffin (e.g. REF 3297 Medical white oil).

### Blunt scissors

Each cutting instrument suffers under abrasion. Scissors which have been properly cleaned should not come into contact with other instruments, especially in the "Ultrasonic bath" or in the "Thermodisinfector", to protect the delicate blades.

We recommend cleaning in washtrays which ensure a fixing of the scissors to avoid contact damage.

**MAKE SURE TO USE YOUR VALUABLE SCISSORS ONLY FOR THEIR INTENDED APPLICATION!**

**FOR MORE INFORMATION ABOUT MAINTAINING AND REPROCESSING OF YOUR INSTRUMENTS, PLEASE VISIT THE HOMEPAGE OF THE WORKING GROUP FOR INSTRUMENT PROCESSING (AKI) UNDER <http://www.a-k-i.org>**

### Recommended care products:

#### 3297

Paraffinum Perliquidum, 50 ml  
(Paraffin Oil / Medical Whiteoil)

Medical Whiteoil for maintenance of instrument hinges and for use as sharpening oil.

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A. Schweickhardt GmbH & Co. KG · Tel.: +49 7464 98910-700 · Fax: -719 · contact@schwert.com · www.schwert.com