

Eliminate cross contamination with Riskontrol

The use of disposable products in the delivery of healthcare is interpreted as assurance to patients that the dental practice is providing a high level of asepsis and cleanliness for the protection of their wellbeing.

It is recognised that the multifunction air/water triplex syringe is used more than any other instrument in the dental surgery. Yet, it has been shown to be virtually impossible to attain full cleaning and sterility of the lumen, which means that it is a high-risk source of cross-contamination. When the syringe is placed intraorally it mixes with oral contaminants such as blood, saliva, bacteria and viruses via direct and aerosol contact: both the inside and outside of the triplex syringe become tainted.



100%
Hygiene
Ergonomy
Dry Air
Ecology
Universal

The multifunction syringe is the instrument most widely used by dentists and constitutes a permanent factor of cross and direct contamination.
A fixed element of your unit, the syringe cannot be sterilized in the autoclave or by any other means, unlike turbines, contra-angles.
Without a disposable tip, your air/water syringe is the weak point of the aseptic chain of your dental surgery unit. The disposable syringe tip also avoids microbial migration in the unit's water and air circuits and therefore ensures compliance with the strictest hygiene standards.

Satelec Pierre Roland of France has introduced a new disposable triplex syringe tip called Riskontrol that eliminates cross-contamination from the most overlooked source in the dental practice.

The Riskontrol Disposable Triplex Syringe Tip was developed from extensive research and development. The tips are used worldwide and retain a multifunctional performance by delivering a strong aerosol spray, or an independent air or water stream. Their unique design incorporates two independent tubes that guide air and water to the extremities of the tip providing totally dry air to the mouth. This is essential for any cementation procedures or operations that require a totally dry field.

Any triplex syringe unit can be converted with the availability of more than 70 different anodised screw or clip-on adaptors. All adaptors come with 2 protective tips and (depending on the type) a replacement seal to ensure optimum syringe performance.

Riskontrol tips are easily fitted by pushing them on to the adapter and with removal as simple as lining up the marks and pulling the tip off. You only have to install the adapters once so no repeated adjustments are necessary.

Riskontrol tips are flexible and bendable. Each tip is made from non-shape memory plastic derived from 100% recyclable food-quality plastic, which makes it an environmentally safe alternative for a disposable product. The bulk 250 quantity dispenser box also keeps them economical.

Riskontrol tips are exclusively distributed by Gunz Dental. Call 1800-025-300 or 0800-301-010 in New Zealand to place an order or request additional information.



ADX...
everything for
the professional
clinical dental
team.

Australian Dental Expo
Melbourne Exhibition Centre
26-28 May 2006
In conjunction with
the ADAVB Convention



FREE registration online now at

www.adia.org.au