

REGISTRATION FORM

Gunz Account No: _____

Name: _____

Name/s to appear on Certificate/s of Attendance:

Address: _____

Email: _____

Attendees will receive confirmation of booking via email.

Tel: _____ Fax: _____

Payment: Cheque (payable to Gunz Dental Pty Ltd)

Credit Card

No: _____

Exp: _____

CCV _____ (last 3 digits on the back of your credit card)

Signature _____

TO REGISTER

FAX MAIL fax form to (FREEFAX) 1800 025 202
detach and mail to:

Gunz Dental,
Locked Bag 5000,
Alexandria NSW 2015

PHONE Call Alice at Gunz Dental on (02) 9935 6677

COURSE DETAILS

Dates & Venue

QLD Course: Wednesday, 29th July, 2009
Gunz Brisbane Showroom
Unit 5, 106 Fison Ave West,
Eagle Farm 4009.

Lecturer: Julie Carranza, RDH, BScDH.

VIC Course: Wednesday, 23rd September, 2009
Gunz Melbourne Showroom
Ground Floor, 192B Burwood Road,
Hawthorn 3122.

Lecturer: Anne Di Paolo, ADOH (DENTAL HYGIENE).

Registration: 6.00pm Welcome & Refreshments

Workshop Time: 6:30pm – 9:30pm

Cost \$30.00 (incl GST)

Each workshop limited to 25 participants.
Please advise of any special dietary requirements and we will
attempt to address these.



TALK TO US



presents

Hygienist to Hygienist

Mastering Instrument Sharpening

'The 1,000k Tune Up!'

with

Julie Carranza RDH, BScDH

& Anne Di Paolo ADOH (DENTAL HYGIENE)



HANDS-ON WORKSHOP INVITATION

Brisbane • Melbourne

Ph: 1800 025 300 Fax: 1800 025 202 www.gunz.com.au



JULIE CARRANZA RDH, BScDH.

Julie Carranza is an American Dental Hygienist working and living in Queensland. She has been a clinician in a variety of settings including both general dentistry and Periodontics since 1991.

Julie is the President of The Dental Hygienists' Association Inc. Queensland branch. She writes a regular column, *Beyond The Pond*, in the monthly

e-newsletter and works full time at a large practice in Caloundra. She has facilitated a dental study group on the Sunshine Coast, taught oral health to nursing staff and has presented at a variety of events.



ANNE DI PAOLO ADOH (DENTAL HYGIENE).

Anne Di Paolo is the President of Dental Hygienists' Association Inc. Victorian Branch. After graduating as a senior Dental Assistant, Anne joined the ADAVB Dental Nurses' Training Committee and later became involved with RMIT lecturing and mentoring Dental Assistants where she also contributed to the publication of training manuals.

Anne has worked in both general and specialist practices as a Dental Hygienist and has gained substantial experience in all aspects of dentistry. She is also involved with the ADAVB on the Dental Health Education Committee and with the University of Melbourne School Advisory Committee. Anne currently works in general practice.

Hygienist to Hygienist

Mastering Instrument Sharpening

'The 1,000k Tune Up!'

COURSE OUTLINE

- 1) Understanding instrument care and management.
- 2) Review of instrument types, blades and handles.
- 3) Introduction to advanced instrumentation.
- 4) Ergonomic considerations.
- 5) Understanding of and proficiency in instrument sharpening.

Participants will be sharpening the instruments they bring to the workshop. Please bring sickles, curettes and graceys.

LEARNING OBJECTIVES

- 1) Explanation of hand instrument manufacturing, infection control, cassettes.
- 2) Explanation and demonstration of hand instrument variations including: rigid gracey, standard gracey, pocket scalers, mini pocket scalers, files, hoes, furcation scalers, implant scalers. Explanation and demonstration of handle types. Explanation and demonstration of calculus detector (CF 11/12).
- 3) Ergonomic considerations including loupes, fulcrums, strokes, utilising assistants.
- 4) Considerable proportion of the lecture will cover hands on practice with instrument sharpening using both PerioStar and D.I.S.C. Participants will be teamed in small groups with individual opportunity and tutorial. Participants will be sharpening the instruments they bring to the workshop.

Both courses are approved for Continuing Professional Development Points by Dental Practice Board of Victoria.



The Periostar® 3000 from Kerr Hawe efficiently and automatically sharpens the complete range of periodontal instruments. With simple operation and standardised sharpening at a constant angle, Periostar provides consistent results for patient comfort and staff efficiency.



The fully autoclavable D.I.S.C. instrument sharpener from Premier quickly sharpens scalers and curettes correctly without damaging the cutting surface. Ideal for routine sharpening, as well as for chariside touch-ups, the D.I.S.C. sharpener is an easy-to-use, efficient way to sharpen instruments for optimum results.