

TOOTH WHITENING

Offering tooth whitening in your practice is a great way of generating appointments and referrals as well as increasing your profitability.

As an increasingly popular procedure, it is important that you have a fully integrated product offering for your patients. Our tooth whitening range contains a full suite of products, ranging from initial tray fabrication through to in-surgery and take home kits. We also offer a range of in-surgery marketing material to help you promote tooth whitening to your patients.

For all your Tooth Whitening needs Talk To Us.



SECTION 16

TOOTH WHITENING

CONTENTS

TRAY FABRICATION	309
TOOTH WHITENING ACCESSORY TRAY CASE	309
TOOTH WHITENING IN SURGERY	310
ENAMEL MICROABRASION	311
TOOTH WHITENING TAKE HOME	311-313
OPAESCENCE® MARKETING MATERIAL	314

TRAY FABRICATION



Ultradent LC Block-Out Resin is applied to provide reservoir space in bleaching tray.

Instant hard, strong, no-mix material for blocking out undercuts on dies and filling in voids.



LC Block-Out Resin
ULTRADENT

A low viscosity, blue pigmented light cured resin designed to provide reservoir space for bleaching trays. Attach black mini tip onto LC Block Out syringe for rapid and efficient delivery. Cure with curing light for 20 seconds.

Kit - 4 x 1.2ml syringes & 20 black mini tips

ULT0240 

Refill 4 x 1.2ml syringes

ULT0241 

Econo Kit 20 x 1.2ml syringes & 20 black mini tips

ULT0242 



Tray Magic Soft Tray Trimmer
PREMIER

Designed especially for soft tray material including soft denture reline, athletic mouth guards and snore guards. The curved, double edged, heated tip is perfect for cutting scalloped trays. One year warranty. 220V.

Tray Magic

PRE8207 



Sof-tray® Classic Sheet
ULTRADENT

Tray fabrication material. Use 0.035" sheets for most trays and the 0.060 " sheets for patients who have severe bruxism.

Sof-Tray Sheets .035" (0.889mm) Pack of 25

ULT0226 

Sof-Tray Sheets .060" (1.52mm) Pack of 20

ULT0227 



Ultra-Trim Scissors
ULTRADENT

Fine serrated stainless steel scissors, used to precisely shape and trim the border around the interdental papilla when "scaloping" trays.

Ultra Trim Scissors

ULT0605 

TOOTH WHITENING ACCESSORY



Opalescence® Tray Case
ULTRADENT

Protects bleaching tray when not in use and fits easily into pocket or purse

Opalescence pocket tray cases pack

ULT0707/20 

TOOTH WHITENING IN SURGERY



Opalescence® Boost PF
ULTRADENT

Opalescence Boost PF now has potassium nitrate and fluoride to help combat sensitivity. It is a chemically activated whitening that offers the convenience of syringe delivery, plus the bleaching power of activated 38% hydrogen peroxide. Total application time is 45 minutes per visit. Mix to activate the product just prior to application. Activated gel is then applied to teeth for whitening - without a light!

Opal Boost PF Intro Kit Contents:

- 4 x 1.2ml Boost & activator syringes
- 20 x Delivery tips
- 2 x 1.2ml OpalDam syringes
- 2x Surgi Tips
- 2 x IsoBlocks
- 2x Luer Vacuum Adapters

Opalescence Boost PF Intro Kit

ULT5330 

Opal Boost Patient Kit Contents:

- 1 x 1.2ml Boost & activator syringes
- 1 x IsoBlock
- 1 x 1.2ml OpalDam syringe
- 10 x Delivery tips

Opalescence Boost PF Patient Kit

ULT5332 



Opalescence® Quick™ PF 45%
ULTRADENT

Highly viscous, sticky, 45% Carbamide peroxide whitening gel with potassium nitrate and fluoride. This will offer an excellent office supervised, tray delivered bleaching result. This is ideal for jump starts and touch ups without the need of gingival protection. Treatment times may vary from 30 to 60 minutes daily.

Opalescence Quick™ PF 45% Mini Kit contents:

- 4 x 1.2ml Opalescence Quick syringes,
- 4 x delivery tips

Opalescence® Quick™ PF 45% Mini Kit

ULT5346 

Opalescence Quick™ PF 45% Econo Kit contents:

- 20 x 1.2ml Opalescence Quick syringes
- 20 x delivery tips

Opalescence® Quick™ PF 45% Econo Kit

ULT5348 



Opalescence® Endo
ULTRADENT

Pre-mixed, single-syringe delivered 35% hydrogen peroxide gel. Opalescence® Endo is designed for “walking bleach” technique. The whitening gel is activated at body temperature and remains active for 3-5 days. It is perfect for whitening endodontically treated tooth.

Opalescence Refill contents:

- 2 x 1.2ml pre-filled syringes

Opalescence® Endo Refill

ULT1323 



OpalDam®
ULTRADENT



OpalDam is a passively adhesive (sealing), methacrylate-based, light-reflecting resin barrier for isolating tissue adjacent to teeth being whitened. For single tooth bleaching, OpalDam protects adjacent teeth.

OpalDam® Kit contents:

- 4 x 1.2ml OpalDam syringes
- 10 x Micro 20g tips
- 10 x Black Mini tips

OpalDam® Kit

ULT0324 



IsoBlock®
ULTRADENT





Offers a bilateral support with tongue restraint and eases TMJ strain during in-office whitening procedure. Designed to be used with Opalescence® X-tra® Boost.

IsoBlock® Pack of 10

ULT0331 

ENAMEL MICROABRASION



1. Opalustre's chemical and mechanical abrasion produces a natural appearing surface.

2. Silicone carbide microparticles contained in Opalustre.

Opalustre®
ULTRADENT

A viscous, water-soluble, chemical and mechanical abrasion slurry containing hydrochloric acid and silicon carbide micro particles. Use to remove unsightly enamel decalcification defects that are less than 0.2mm in depth. This is ideal for superficial white and brown "demineralisation" due to enamel mottling from fluorosis.

Opalustre Kit contents:

- 4 - 1.2ml Opalustre syringes,
- 20 x White Mac tips,
- 20 x OpalCups (10 bristle & 10 finishing)

Opalustre Kit

ULT0554



Prema
PREMIER

Prema microabrasion is a patented acid abrasion formulation used to correct superficial white, brown or multicoloured "demineralisation" enamel defects.

Prema Kit: 10 cups, 2 mandrels, 10 mandrel tips, 1 finger ring.

PRE8009

TOOTH WHITENING TAKE HOME



REV!
PREMIER

REV! whitening treatment is a revolutionary process that delivers optimal whitening in only 15 minutes, once a day. 14% Hydrogen Peroxide.

Whitening Refresher Packs x 6
PRE0141

Whitening Patient Packs x 6
PRE0142



Bravo
PREMIER

Shorter wear time only 7 hours is required to complete the 2 week take home treatment. No sensitivity patented water-based 9% hydrogen peroxide gel formulation hydrates teeth as they whiten to avoid sensitivity. Cosmetic patient pouch reflects well on your practice, and organises educational material allowing you to send home other practice building items.

Patient Pack 4x3ml Syringe
PRE0092

Value Pack 24 x 3 ml Syringe, 6 Travel Pouches and 6 patient guides.
PRE0093

Touch-Up Pack 1 x 3ml Syringe, dispensing tips.
PRE0091

TOOTH WHITENING TAKE HOME



Opalescence PF
ULTRADENT


It's hard to beat science! Exciting findings of the Opalescence PF product line presented at the American Association of Dental Research (AADR) in 2006 confirm that whitening with Opalescence PF will provide ANTI-Caries benefits, Increase Enamel Microhardness, Minimise sensitivity and Improve Enamel Health.

Opalescence PF is highly recommended for whitening discoloured teeth prior to placement of composite, veneers, and/or crowns. Opalescence PF is effective in removing some or all of the internal tooth discolourations due to congenital, systemic, pharmacologic, traumatic, etc., factors and aging. Great success with fluorosis and even tetracycline staining.

The Opalescence PF sticky formula does not leach from the tray like other bleaching gels. Patients experience results quickly, which motivates them to pursue the treatment. Vital teeth whitened with Opalescence remain lighter for years. Opalescence PF is available in 4 concentrations, 4 formulations, 3 flavours and 2 kit configurations to meet all your whitening needs.



Sticky, viscous formula does not leach from the tray like other bleaching gels.



Indiana University: Small quantities of gel were removed and analysed at various intervals.

Doctor's Kits
ULTRADENT

An impressive presentation to patients and ideal for surgeries who purchase tray material in bulk or have a lab fabricate their trays

Doctors Kit contents:

- 8 x 1.2 ml prefilled syringes
- 1 x Tray case
- 1 x Cosmetic Bag
- 1.2 ml LC Block out
- 2 x 0.035" 5"x 5" Sof-tray sheets
- 1 x Shade Guide

10% Regular PF
ULT5381

10% Mint PF
ULT5379

10% Melon PF
ULT5380

15% Regular PF
ULT5384

15% Mint PF
ULT5382

20% Regular PF
ULT5387

20% Mint PF
ULT5385

20% Melon PF
ULT5386

35% Regular PF
ULT5390

35% Mint PF
ULT5388



Refill Kits

Suitable for surgeries that purchase materials in bulk and assemble their own kits.

Contents: 40 x 1.2ml prefilled syringes

10% Regular PF
ULT5396

10% Mint PF
ULT5394

15% Regular PF
ULT5399

15% Mint PF
ULT5397

20% Regular PF
ULT5402

20% Mint PF
ULT5400

20% Melon PF
ULT5401

35% Regular PF
ULT5405

35% Mint PF
ULT5403



TOOTH WHITENING TAKE HOME



Ultra EZ® Desensitisers
ULTRADENT

UltraEZ® is ideal to use to quickly eliminate sensitivity from toothbrush abrasion, thermal and chemical changes. As well as tooth bleaching and root exposure. It is available in syringe to use with custom trays and disposable pre-loaded trays. UltraEZ® contains Potassium Nitrate and Fluoride to address both short and long term sensitivity.

UltraEZ Syringe Kit

ULT1008 

UltraEZ Tray Delivered Upper

ULT7200 

UltraEZ Tray Delivered Lower

ULT7201 



Opalescence Trèshwhite® Supreme
Ultradent



New Opalescence Trèshwhite Supreme is the first pre-loaded, convenient tray delivered over-the-counter take home whitening system in the market. Opalescence Trèshwhite Supreme contains 10% hydrogen peroxide for faster and superior whitening results. It also contains potassium nitrate and fluoride. The convenience of pre-loaded trays of Opalescence Trèshwhite Supreme means no lab and tray fabrication is required. This product is also ideal for patients where impression taking is a challenge. The sleek design of the trays provides ease of use and comfort to the patient. The external and internal bleaching trays can easily be positioned and adapted to fit any smile. Opalescence Trèshwhite Supreme offers 3 delicious, improved flavours: mint, melon and peach. The optimum gel quantity of this product allows for rapid, easy cleanup following bleaching. The new high-tech packaging of Trèshwhite provides superior stability and shelf life. Opalescence Trèshwhite Supreme PF outstanding unique formula will give patients fresher and whiter teeth while rebuilding enamel and preventing caries. All Patient Kits Contain 10 upper and 10 lower pre-filled trays.

Value Pack Mint - Pack of 6 Kits

ULT5706-Value Pack



Value Pack Melon - Pack of 6 Kits

ULT5707-Value Pack



Patient Kit Mint

ULT5706



Patient Kit Melon

ULT5707



OPALESCENCE® MARKETING MATERIAL*

A full range of in surgery marketing material is available. Items shown are for illustration purposes. Contact your Gunz representative for a sample.



Posters

Made to fit frames, they make a great statement when displayed in surgeries and waiting rooms.



Ceiling Tiles

Place these on the ceiling above dental chairs. Designed to excite patient interest in whitening procedures.



Waiting Room Brochures

Offer a quick and informative Q & A format on Tooth Whitening as well as before and after photographs.



Recall Cards

Send as post card appointment reminders to your patients.

*** Actual artwork may vary subject to availability.**