

SECTION 7

Ensuring that your surgery maintains the best infection control standards is an ongoing and costly exercise. Gunz Dental stocks a complete range of products to satisfy all requirements for both you and your patients.

For all your Infection Control needs Talk to Us.

INFECTION CONTROL



SECTION 7

INFECTION CONTROL

CONTENTS

BARRIER PRODUCTS	125-126
AIR WATER SYRINGE TIPS	126
CLEANERS	127
DETERGENTS	127-128
DISINFECTANTS	128
MASKS	128-129
GLOVES	129
HANDWASHES	130
EYEWEAR	131
GAUZE	131
MONITORING STRIPS	131-132
STERILISATION POUCHES	132
SCALPEL BLADE REMOVER	132

BARRIER PRODUCTS

A range of barrier products designed to cover all essential items used in your day to day procedures. They provide adequate infection control protection and peace of mind for your patient.



Tray Sleeve
SDS KERR

10.5" x 14" Pack 500
KER3000-B1

11.5" x 16" Pack 500
KER3200-E



Tube Sleeve
SDS KERR

2" Diameter 1200 foot roll
KER3720



Syringe Cover
ULTRADENT

A thin, textured, polypropylene cover designed to custom-fit Ultradent 1.2ml syringes with tips attached.

Syringe Covers Pack 250
ULT0249



Head Rest Covers
SDS KERR

Designed to preserve the aesthetic appearance of the dental chair and provide protection against body oils and perspiration.

11.5" x 9.5" Pack 250
KER3511

14" x 9.5" Pack 250
KER3514



PROGRESS

White 2ply 25x25cm Carton 500
MDC3014

White 2ply 25x33cm Carton 500
MDC3017



Light Sleeve
SDS KERR

Fits T Light handles 5.75" x 4" Pack 500
KER3600



Syringe Sleeve
SDS KERR

Disposable sleeves made to cover the 3 way syringe, HVE and saliva ejector valves.

Air Water Syringe 2.5" x 10" Pack 500
KER3760-1

With Adhesive 2.5" x 11.5" Pack 500
KER3763-WA1

With Opening 2.5" x 10" Pack 500
KER3767-W01



Handpiece Sleeve
SDS KERR

Straight Pack 500
KER4000

Contra Angle Pack 500
KER4050



Chair Sleeve
SDS KERR

27.5" x 24" Pack 225
KER3800

CLEANERS



Vac Attak™
PREMIER

High performance evacuation system cleaner, especially effective in removing prophy paste and fluoride gels, the two most common causes of blockages.

800 gram jar with scoop dispenser

PRE11101 



Puli-Jet
CATTANI

A disinfectant solution especially formulated for the cleaning and deodorising of dental aspirator lines. Its active ingredients have been specifically selected to ensure all plastic and rubber components remain undamaged in the cleaning process.

1 Litre

CLK0065 

5 Litres

CLK0065/5 

Puli-Jet Pump for 5 Litre

CLKPJP/5 

10 Litre

CLK0065/10 



Isopropyl Alcohol
DENTALIFE

A general purpose alcohol that provides antibacterial, cleaning and solvent properties in the surgery. Comprised of 70% Isopropyl Alcohol and 30% purified water it is suitable for cleaning, surface and instrument preparation.

Clear 70%

5 Litre

DL0901 

Blue 70%

5 Litre

DL0911 



Vacusan
PROGRESS

A non toxic, non corrosive enzyme-based cleaner, designed specifically to clean, sanitise and deodorise internal components of all dental surgery suction equipment.

5 litres

PRO900 



Microvac
WHITELEY

High powered cleaning ability in a non-foaming formulation that eliminates foul odours.

1kg

WHI130194 

DETERGENTS



Mediclean
WHITELEY

A non-toxic, hospital grade neutral detergent that provides a streak free finish on all surfaces. It is a clear liquid with a light lemon fragrance. Coupled with Viraclean, Mediclean makes for simplification of the cleaning and disinfection process.

500ml

WHI130204 



Sonidet
WHITELEY

The most widely used neutral instrument detergent in Australia. Ideal for use in ultra-sonic machines and the manual cleaning of instruments, Sonident is an integral part of the sterilisation process. It is odourless, non staining and is a safe detergent for the removal of all surgical soils.

5 Litre

WHI190377 



Ultrasoak
PROGRESS

A specially formulated dual-purpose neutral detergent concentrate used as an immersion solution to help remove materials from instruments.

5 Litres

PRO905 

DETERGENTS



Clinicare
DENTALIFE

Neutral detergent products formulated to address tough clinical cleaning situations.

Dispenser Pack - Contains 50 x 30ml Sachets

DL2605 

Canister 220 Towelettes x 12

DL2610 

Towelette Towelette Tub 220 wipes

DL2610Z 

Refill Pack 220 wipes

DL2610RZ 

DISINFECTANTS



Aidal Plus
WHITELEY

A new and improved instrument disinfectant designed with 2% Glutaraldehyde for use in a surgical practice. It is deodorised.

5 Litre

WHI010082 

DISINFECTANTS



Bactol Blue
WHITELEY

A ready to use hospital grade disinfectant specifically designed for the rapid killing of pathogenic microorganisms.

5 Litre

WHI020105 



Viraclean
WHITELEY

The most popular hospital grade surface disinfectant used in health care in Australia. A unique patented formulation suitable for disinfection of surfaces soiled by blood or body fluids.

5 Litre

WHI210556 

500ml Squeeze bottle

WHI210574 

MASKS



MEDICOM

Earloop masks are very soft, don't pull on the ear, are non fluffy, odour free and do not contain fiberglass.

Tie on Blue Box 50

MDC2000 

Earloop Blue Box 50

MDC2015 



SMITH & NEPHEW

Duckbill Surgical Tie on Box 50

SN1020K 

Earloop Splash Resistant Box 50

SN103414 

High Filtration Splash Resistant Tie on Box 50

SN211666 



3M ESPE

Features slip-on convenience and provides fluid resistance and a bacterial filtration efficiency of 99% at a mean particle size of 4.0 microns.

Fluid Resistant Box 50

MMM1942FB 

MASKS



3M ESPE

Known as the industry standard this comfortable mask provides proven fluid resistance.

MMM1818 Tie-On Box 50
3M ESPE



Moulded
3M ESPE

Top of the line mask with exceptional fluid resistance. The unique triple layer design includes a fluid resistant outer layer, a micro fibre middle layer that traps micro organisms and an absorbant inner layer which absorbs moisture to provide comfort.

Triple Layer ear loop Box 50
MMM2042F

Ear Loop Box 50
MMM2000F



Fluidshield Fog Free Procedure Mask
KIMBERLY CLARK

Tecnol Fluidshield fluid resistant face mask meets the industry's highest standard for facial protection against the penetration of blood and body fluids. PFE 99% & BFE 99%

Pleat-style with ear loops Box 40
KIM47107

GLOVES



Green - Large
Blue - Medium
Pink - Small
Purple - X-Small

PROGRESS

Made of natural latex with a low protein content, minimising the occurrence of latex allergies. Available in both lightly powdered and powder free formulations.

Box- Contains 100 Gloves

Low Powder	X-small	PROLPXSZ	
	Small	PROLPSZ	
	Medium	PROLPMZ	
	Large	PROLPLZ	
Powder Free	X-small	PROPFXSZ	
	Small	PROPF SZ	
	Medium	PROPFMZ	
	Large	PROPFLZ	

Carton- Contains 10 Boxes

Low Powder	X-small	PROLPXS	
	Small	PROLPS	
	Medium	PROLPM	
	Large	PROLPL	
Powder Free	X-small	PROPFXS	
	Small	PROPFS	
	Medium	PROPFM	
	Large	PROPFL	

HANDWASHES

Avagard Hand and Body Washes

3M ESPE

All antiseptic washes are broad spectrum bactericides for use as procedural handwashes. The general wash is compatible with and does not reduce the efficiency of the Avagard antiseptics. All are skin compatible with a soap free high level moisturising base.



General 500ml
MMM9230-P



Antiseptic with 2% Gluconate 5 litre
MMM9232-B



Avagard Surgical Hand Scrubs

3M ESPE

Antiseptic Surgical handscrubs with broad spectrum bactericide. Exceptional skin compatibility.

With 0.75% Povidone -iodine 125ml
MMM9240-H

With 0.75% Povidone -iodine 500ml
MMM9240-P

With 0.75% Povidone -iodine 1.5 Litre
MMM9240-D



With 4% Chlorhexidine Gluconate 500ml Bottle
MMM9241-P



Avagard Moisturising lotions

3M ESPE

High level moisturising emollient barrier lotions. Leaves protective film restricting moisture loss. Non-greasy, non-slippery and contains no silicone. Compatible with latex gloves and Avagard antiseptics.

General 125ml
MMM9260-H

General 500ml
MMM9260-P

Antiseptic with Triclosan 125ml
MMM9261-H

Antiseptic with Triclosan 500ml
MMM9261-P

Antiseptic with Triclosan 1.5 Litres
MMM9261-D



Avagard Handrubs

3M ESPE

Waterless, alcohol-based handrubs. Topical broad spectrum bactericide. Residual and cumulative activity. Fast drying with good skin compatibility.

With 0.5% Chlorhexidine Gluconate 125ml
MMM9250-H

With 0.5% Chlorhexidine Gluconate 500ml
MMM9250-P



Microshield Handwash

JOHNSON & JOHNSON

Microshield Handwash removes debris and transient microorganisms. Moistures, with emollients to minimise the potential drying and skin irritation effects associated with frequent handwashing.

500ml Box of 12
JJM61392

1.5L Box of 6
JJM61224

5L Box of 2
JJM60238

EYEWEAR



Safety Glasses
ULTRADENT

High quality, light weight protective eyewear with adjustable armpiece. Moulded polycarbonate to ensure high optical properties and scratch resistance. Stylish wrap around design for side protection.

Smoke Frame with Smoke Lens

ULT0914 

Clear Frame with Clear Lens

ULT0501 

GAUZE



Multisorb
SMITH & NEPHEW

Non woven 7.5 x 7.5cm Pack 100
SN2307 



Topper 8
JOHNSON & JOHNSON

High absorbency without linting, soft & strong with minimal allergenic potential and provides a textured surface to facilitate debris removal.

Non-sterile:

5 x 5cm Pack of 100
JJMM50405 

7.5 x 7.5cm Pack of 200
JJMM53407N 

10 x 10cm Pack of 100
JJMM50410 

Sterile:

7.5 x 7.5cm Pack of 100
JJMTS8075 

5 x 5cm Pack of 25
JJMTS8055 



Safe Gauze
MEDICOM

Non woven, non-sterile 4 ply sponges.
Highly absorbant and non-linting.

7.5 x 7.5 cm Pack 200
MDC2103-CH 

MONITORING STRIPS



Comply™
3M ESPE

These indicators are placed inside a pouch or instrument cassette and indicate whether the sterilant has penetrated the pack. The colour changes from off white to black when fully processed. It does not verify sterility but allows detection of certain procedural errors and equipment malfunctions.


Comply Chemical Indicator Strips Refill

MMM1250 



Attest™
3M ESPE

The indicator vial is placed in or on the autoclave pouch or instrument cassette and measures the microbial killing power of a sterilisation process. The Blue indicators are used in benchtop sterilisers where the come up time to temperature is equal to or less than 1 minute. The Brown indicators are for use in sterilisers where the come up time is greater than 1 minute. In sterilisers where both wrapped and unwrapped items are being processed, the Brown indicator is the correct one to use.

Attest Biological Indicator Blue Pack 25
MMM1261P 


Attest Biological Indicator Brown Pack 25
MMM1262P 

MONITORING STRIPS



Bowie Dick
3M ESPE

A product to assess air removal and steam penetration in pre vacuum autoclaves. It does not prove sterility but shows whether the vacuum system of an autoclave is operating correctly. It is recommended that all class B autoclaves be tested on a daily basis as part of normal Quality Assurance measures.

3M Bowie Dick Test Carton of 20
MMM1300 

3M Bowie Dick Test x 1
MMM1300Z 

STERILISATION POUCHES

85 x 159mm Box 200
MDC68005 

89 x 254mm Box 200
MDC68010 

57 x 130mm Box 200
MDC68015 

134 x 279mm Box 100
MDC68025 

254 x 381mm Box 100
MDC68035 


SCALPEL BLADE REMOVER




Qlicksmart
STENNING

Designed for single-handed operation, this device combines the process of removal, containment and disposal of scalpel blades in one simple procedure. Once removed the blade has no further human contact. It is designed for use with most # 3, #4, # 7 and #9 scalpel blades. It comes with a progressive counter and will remove a maximum of 100 blades and includes an automatic shut off mechanism to prevent over filling.

A universal wall bracket is available separately that can be used to mount the device on a wall, bench or trolley.

Qlicksmart Flask Scalpel Blade Remover
QS-FS-1 

Qlicksmart Universal Bracket
QS-BR-1 

STERILISATION POUCHES



Medicom

See the difference with built in colour change indicators along with puncture resistant, blue tinted see-through plastic. Extra strength provided with heavy medical grade paper.

Compatible with steam, chemical vapour and gas methods.

Brands you can trust

- ★ 3M ESPE
- ★ ACTEON GROUP
- ★ CATTANI
- ★ DENTALIFE
- ★ JOHNSON & JOHNSON
- ★ KIMBERLY CLARK
- ★ MEDICOM
- ★ PREMIER
- ★ PROGRESS
- ★ SDS KERR
- ★ SMITH & NEPHEW
- ★ STENNING
- ★ ULTRADENT
- ★ WHITELEY

For all your infection control needs

TALK TO US