

PREVENTATIVE

The aim of all dental practitioners is to help patients keep their teeth and gums healthy every day. Our preventative product range includes items for cleaning, flossing and polishing your patients' teeth as well as items for ensuring good oral hygiene practices.

For all your Preventative needs Talk To Us.



SECTION 10

PREVENTATIVE

CONTENTS

DISCLOSING SOLUTION	185
MOUTH RINSE	185
FLOSS	185
PIT AND FISSURE	185-186
PROPHY BRUSHES	186
PROPHY CUPS	186-187
PROPHY PASTE	188
SCALER	188
DENTURE CLEANER	188
TOOTHBRUSHES	188
MOUTH MOISTENER	189
INTERDENTAL BRUSH	189

DISCLOSING SOLUTION



Red-Cote
PREVIDENT

Highlights harmful bacterial plaque on tooth surfaces and shows patient areas where more brushing and flossing are needed. Non toxic cherry flavour.

Bottle 6cc

PVD806 

Tablets Pack 1000 Cello wrapped tablets

PVD800 

MOUTH RINSE



Swirl
DENTALIFE

This mouth rinse concentrate cleans, clears and refreshes the mouth after every rinse. The economical formulation dispenses up to 1300 applications in each bottle.

Wildberry Bottle 200ml

DL0811 

Freshmint Bottle 200ml

DL0801 

FLOSS



Dentotape
JOHNSON & JOHNSON

5 yards, Pack 144

DPA35312 

100 yards, Pack 6

DPA35316 

20 metres, Pack 6

DPA35701 

Waxed
JOHNSON & JOHNSON

200 metres, Pack 4

DPA35700 

Mint 50 metres, Pack 6

DPA35702 

50 metres Pack 6

DPA35704 



Eez-Thru Floss Threaders
BUTLER

Versatile nylon loop helps patients thread floss through spaces where most other products cannot reach. Suitable for use under crowns and around orthodontic and implant appliances

Pack 25

PVD840 

Pack 5

PVD840P 



Flossmate Handle
BUTLER

Longer curved fingers allow access to posterior teeth and make it easy to place floss between teeth without interference from the handle.

Flossmate Floss Handle

PVD845 

PIT & FISSURE



UltraSeal XT® plus™
ULTRADENT

The superhero in the battle against decay. UltraSeal XT plus pit and fissure sealant provides exciting improvements in sealant technology. The light-cured, 58% filled resin makes UltraSeal XT plus a stronger, more wear-resistant composite sealant. UltraSeal XT plus has less polymerization shrinkage than competitive products. The spiral brush action of the Inspiral Brush tip causes shear thinning of the filled, thixotropic UltraSeal XT plus, therefore decreasing its viscosity as it is placed. Resin firms when shear thinning ceases (placement complete) and helps prevent resin from running before it can be light cured. UltraSeal XT plus is significantly superior to other pit and fissure sealant materials.

UltraSeal XT plus kit contents:

- 1 x 1.2ml UltraSeal XT plus syringe
- 1 x 1.2ml UltraEtch syringe
- 2 x 1.2ml PrimaDry syringes
- 20 x Inspiral Brush tips
- 20 x Blue Micro Tip

UltraSeal XT plus Refill kit contents:

- 4 x 1.2ml UltraSeal XT plus syringes

PIT & FISSURE

UltraSeal XT plus (Opaque shade only) Econo Refill contents:

20 x 1.2ml UltraSeal XT plus syringes

UltraSeal XT plus Kit Shade Clear

ULT0563 

UltraSeal XT plus Kit Shade Opaque White

ULT0725 

UltraSeal XT plus Kit Shade A2

ULT0733 

UltraSeal XT plus Refill Shade Clear

ULT0565 

UltraSeal XT plus Refill Shade Opaque White

ULT0726 

UltraSeal XT plus Refill Shade A2

ULT0734 

Opaque White Econo Refill 20 x syringes

ULT0727 

To maximise thixotrophy of Ultraseal XT Plus use Inspiral Brush Tip. See General section page 93.



PrimaDry
ULTRADENT

Use PrimaDry in fissures just prior to resin placement, followed by air drying. PrimaDry enhances resin penetration significantly.

PrimaDry Refill - 4 x 1.2ml PrimaDry Syringes


ULT0716 



Universal Dentin Sealant
ULTRADENT

Nothin' gets past this stuff.

Universal Dentin Sealant is a biocompatible, non-polymerizable, high molecular weight resin in a



BEFORE AFTER

volatile organic solvent. This material is ideal for coating sensitive roots following root planing and minimizing patient post-op discomfort by sealing dentinal tubules. Universal Dentin Sealant forms a virtually impenetrable temporary biocompatible barrier against bacteria, resins, acids, cements, restorations, etc.

Universal Dentin Sealant Refill, pack of 4 x1.2ml syringe

ULT0266 



LuxaFlow Fluorescence
DMG

Ideal for flowable composite Class I's, blocking out undercuts and as a solution for enamel defects. Also effective in creating a seal in the gingival margin interfaces. Delivered from a high-precision syringe tip dispensing system, LuxaFlow Fluorescence provides excellent flowability in a non-slumping formula.

Intro Kit

RG3100 includes 4 Syringes of 1.5 g paste (Shades A1, A2, A3.5, B1) and 20 Luer-Lock-Tips

LuxaFlow Fluorescence Refill A1

RG3100 

LuxaFlow Fluorescence Refill A2

RG3101 

LuxaFlow Fluorescence Refill A3.5

RG3102 

LuxaFlow Fluorescence Refill B1

RG3104 

PROPHY BRUSHES



PROGRESS

Screw-In Junior Cup Pack 100
PRO104 

Screw-In Pack 100
PRO105 

Snap-On Bristle Brushes Pack 100
PRO106 



ICB® Intracoronary Bristle Brush
ULTRADENT

The ICB® brush was designed to clean intracoronally as well as extracoronally before bonding fillings, inlays, onlays and crowns. It is great for cleaning areas with restricted access. It is autoclavable and designed for use with a slow speed handpiece.

Pack 30
ULT1076 
Pack 100
ULT1077 

PROPHY CUPS



PROGRESS

White web design

Screw-in Pack 100
PRO206 

Snap-On Pack 100
PRO207 

PROPHY CUPS



PVD

Screw Type Soft Blue Rib & Web Pack 144
PVD85132

Latch Type Soft Blue Rib & Web Pack 144
PVD85232



PD

Splash proof with contra angle shank

RA White Webbed Pack 12
PRO202

Brush Type RA Pack 12
PRO204



Pro-Cup®
HAWE

The innovative splash free prophy cup with unique lamellar profile provides an economical, splash free application of prophy paste. Thanks to a steady supply of paste, a thin film is built up between the cup and the tooth, reducing heat through friction, providing extra comfort for your patient.

Pro-Cup Standard, Latch-type Pack 30
HAW990

Pro-Cup Hard, Latch-type Pack 30
HAW991

Pro-Cup Screw-type Pack 30
HAW1810

PROPHY PASTE



Optum
DENTALIFE

Is a medium grit, high density, splatter free formula suitable for mild abrasion and stain removal.

Freshmint 200g

DL0201

Banana 200g

DL0202

Bubblegum 200g

DL0203

Tropical 200g

DL0204



Optum F
DENTALIFE

Has all the characteristics of Optum but with 0.22% sodium fluoride.

Freshmint 200g

DL0212

Tropical 200g

DL0214

Wildberry 200g

DL0211



HAWE

Plastic ring pot for single dose application. Fits snugly on finger.

Ring Container Pack 100

HAW985

Prophy Paste Single Dose Box 200

HAW3140



ENAMEL PRO
PREMIER

Creates a lustrous surface!

The only prophy paste with ACP (Amorphous Calcium Phosphate).

It creates ACP when it contacts the patient's teeth and saliva. As ACP forms, it is incorporated within the enamel surface where it remains after rinsing. The available ACP helps prevent future damage by stimulating re-mineralisation of the tooth enamel.

Each Box of Enamel pro contains 200 single-dose cup, packaged in a resealable foil bag and two ComFy-Grip prophy grips to conveniently hold the individual prophy cups.

Mint Fine

PRE7600

Mint Medium

PRE7601

Mint Coarse

PRE7602

Cinnamon Medium

PRE7605

Strawberry Fine

PRE7607

Strawberry Medium

PRE7608

Strawberry Coarse

PRE7609

Grape Medium

PRE7611

SCALER



PREVIDENT

Proflex flexible scaler for calculus removal.

Densco Proflex Scaler 36cm roll

PVD85000



DENTURE CLEANER



Pro Clean
PREMIER

Pro Clean is a professional strength tartar and stain remover for dentures, partials, retainers and night guards. Use in office and dispense for patient home use. A brief 2 hour soak once a week is all that is needed.

Box 25 sachets

PRE11103



TOOTHBRUSHES



BUTLER

End Tuft - Tapered Trim Soft 3 row 7 tufts

PVD308



Classic Soft, 3 Row Professional Pack 6

PVD311



Toothbrush 4 Row Adult Compact UltraSoft

PVD449



TOOTHBRUSHES



Sonicare Electric Toothbrush

PHILIPS



Sonicare's patented technology creates a unique dynamic fluid cleaning action, making it exceptionally effective in hard-to-reach areas. Sonicare delivers a unique oral care experience with results you can see and feel. It is proven to remove more plaque biofilm beyond the bristles and to be 4 times gentler on dentine than oscillating/rotating toothbrushes.

The Sonicare Elite has a 30-second Quadpacer, a 2-minute Smart timer, a deluxe recharge gauge, dual-speed control, a contoured brush head that lasts 6 months and an easy-start function that gently increases brushing power after the first 14 days.

Each pack contains:

- 1 x Toothbrush with 1 head
- 1 x Charger base with easy-wrap cord
- 1 x Brush Head Holder
- 1 x Deluxe Holder
- 2 Year Warranty

Trial Pack for the Dental Professionals
(not for resale)

BREHX7841



Professional Dispensing Pack

BREHX7881



Sonicare Elite replacement head 1/pack

BREHX7001



Sonicare Elite replacement head 2/pack

BREHX7002



Mini Toothbrush Replacement Head Single Pack

30% smaller than the standard brush for focused, attentive and precise cleaning between teeth and in hard-to-reach areas

BREHX7011



MOUTH MOISTENER



Dentacal
NSI

Helps moisten and freshen the mouth while helping fight tooth decay.

Suitable for anyone with a dry mouth whatever the cause, temporary or permanent. Suitable for people with Sjogren's syndrome, chronic xerostomia, damaged salivary glands, or diabetes.

100ml Spray pack

NIUL11790

INTERDENTAL BRUSH

Proxabrush

This system of interdental brushes are designed for easy, convenient removal of plaque from hard to reach tooth surfaces or wide spaces. They are ideal for interproximal spacing, furcations, bridges, extractions as well as orthodontic appliances and implants.



Handle
BUTLER

Double ended handle fits 2 refills at one time for effective cleaning of different size spaces. Its patented lock design allows easy loading and holds refills securely.

Plastic Handle

PVD605




Metal Handle
BUTLER

Textured round handle allows precise control under wet and dry conditions. Screw top closure holds refills securely. Handle is autoclavable.

Metal Single-ended Handle

PVD610



Refill
BUTLER

Nylon filaments and nylon coated wire core provide comfort and will not cause galvanic sensitivity nor scratch implants.

Refill Fine Cylindrical Pack 8
PVD612

Refill Fine Tapered Pack 8
PVD614

Brands you can trust

- ★ BUTLER
- ★ DENTALIFE
- ★ DMG
- ★ HAWE
- ★ JOHNSON & JOHNSON
- ★ NSI
- ★ PD
- ★ PHILIPS
- ★ PREMIER
- ★ PREVIDENT
- ★ PROGRESS
- ★ PVD
- ★ ULTRADENT

For all your preventative needs

TALK TO US