

TOOTH WHITENING

Offering toothwhitening in your practice is a great way of generating appointments and referrals as well as increasing your profitability.

As an increasingly popular procedure, it is important that you have a fully integrated product offering for your patients. Our tooth whitening range contains a full suite of products, ranging from initial tray fabrication through to in-surgery and take home kits. We also offer a range of in-surgery marketing material to help you promote tooth whitening to your patients.

For all your Tooth Whitening needs Talk To Us.



SECTION 14

TOOTH WHITENING

CONTENTS

TRAY FABRICATION	261
TOOTH WHITENING ACCESSORY TRAY CASE	261
TOOTH WHITENING IN SURGERY	261-262
ENAMEL MICROABRASION	262
TOOTH WHITENING TAKE HOME	263-264
OPAESCENCE® MARKETING MATERIAL	265

TRAY FABRICATION

TOOTH WHITENING IN SURGERY



Tray Magic Soft Tray Trimmer
PREMIER

Designed especially for soft tray material including soft denture reline, athletic mouth guards and snore guards. The curved, double edged, heated tip is perfect for cutting scalloped trays. One year warranty. 220V.

Tray Magic
PRE8207



LC Block-Out Resin
ULTRADENT



The utility resin.

A low viscosity, blue pigmented light cured resin designed to provide reservoir space for bleaching trays. Attach black mini tip onto LC Block Out syringe for rapid and efficient delivery. Cure with curing light for 20 seconds.

Kit - 4 x 1.2ml syringes & 20 black mini tips
ULT0240



Refill 4 x 1.2ml syringes
ULT0241



Econo Kit 20 x 1.2ml syringes & 20 black mini tips
ULT0242



FOR YOUR TRAY IMPRESSIONS USE STATUS BLUE ALGINATE - EASY AND REPOURABLE. FOR DETAILS SEE IMPLANT SECTION PAGE 99 OR PROSTHODONTIC PAGE 193.



Sof-tray® Classic Sheet
ULTRADENT

Tray fabrication material. Use 0.035" sheets for most trays and the 0.060 " sheets for patients who have severe bruxism.

Sof-Tray Sheets .035" (0.889mm) Pack of 25
ULT0226



Sof-Tray Sheets .060" (1.52mm) Pack of 20
ULT0227



Ultra-Trim Scissors
ULTRADENT

Not your kindergarten scissors!

Fine serrated stainless steel scissors, used to precisely shape and trim the border around the interdental papilla when "scalloping" trays.

Ultra Trim Scissors
ULT0605



TOOTH WHITENING
ACCESSORY



Opalescence® Tray Case
ULTRADENT

We've made tooth whitening a little safer for the trays.

Protects bleaching tray when not in use and fits easily into pocket or purse

Opalescence pocket tray cases pack
ULT0707/20



Opalescence® Xtra® Boost™
ULTRADENT



No Light, Super Charged, Power White!

Opalescence Xtra Boost offers the convenience of syringe delivery, plus the bleaching power of activated 38% hydrogen peroxide. Total application time is 45 minutes per visit. Mix to activate the product just prior to application. Activated gel is then applied to teeth for whitening - without a light!

Opal Xtra Boost Intro Kit Contents:
4 x 1.2ml Boost & activator syringes
2 x 1.2ml OpalDam syringes
2 x IsoBlocks
20 x Delivery tips

Opalescence® X-tra® Boost Intro Kit
ULT0387

Opal Xtra Boost Mini Kit Contents:
1 x 1.2ml Boost & activator syringes
1 x 1.2ml OpalDam syringe
1 x 1.2ml UltraEZ syringe
10 x Delivery tips

Opalescence® X-tra® Boost Mini Kit
ULT0388



Opalescence® Quick™
ULTRADENT



Whitening from the comfort of your waiting room!

Highly viscous, sticky, 45% Carbamide peroxide whitening gel with potassium nitrate and fluoride. This will offer an excellent office supervised, tray delivered bleaching result. This is ideal for jump starts and touch ups without the need of gingival protection. Treatment times may vary from 30 minutes to 2 hours daily

Opalescence Quick Mini Kit contents:
4 x 1.2ml Opalescence Quick syringes,
4 x delivery tips

Opalescence® Quick™ Mini Kit
ULT5346



TOOTH WHITENING IN SURGERY

Opalescence Quick Econo Kit contents:
20 x 1.2ml Opalescence Quick syringes
20 x delivery tips

Opalescence® Quick™ Econo Kit

ULT0711 



Opalescence® Endo
ULTRADENT

Bringing teeth back from the dead.

Pre-mixed, single-syringe delivered 35% hydrogen peroxide gel. Opalescence® Endo is designed for "walking bleach" technique. The whitening gel is activated at body temperature and remains active for 3-5 days. It is perfect for whitening endodontically treated tooth.


Opalescence Refill contents:

2 x 1.2ml pre-filled syringes

Opalescence® Endo Refill

ULT1323 



OpalDam®
ULTRADENT 

Isolation, protection makes your patient happy.

OpalDam is a passively adhesive (sealing), methacrylate-based, light-reflecting resin barrier for isolating tissue adjacent to teeth being whitened. For single tooth bleaching, OpalDam protects adjacent teeth.

OpalDam® Kit contents:

- 4 x 1.2ml OpalDam syringes
- 10 x Micro 20g tips
- 10 x Black Mini tips

OpalDam® Kit

ULT0324 



IsoBlock®
ULTRADENT

Comfort for you and the patient, who knew?

Offers a bilateral support with tongue restraint and eases TMJ strain during in-office whitening procedure. Designed to be used with Opalescence® X-tra® Boost.

IsoBlock® Pack of 10

ULT0331 

ENAMEL MICROABRASION



Opalustre®
ULTRADENT 

When it's time to break out the heavy equipment!

A viscous, water-soluble, chemical and mechanical abrasion slurry containing hydrochloric acid and silicon carbide micro particles. Use to remove unsightly enamel decalcification defects that are less than 0.2mm in depth. This is ideal for superficial white and brown "demineralisation" due to enamel mottling from fluorosis.

Opalustre Kit contents:

- 4 - 1.2ml Opalustre syringes,
- 20 x White Mac tips,
- 20 x OpalCups (10 bristle & 10 finishing)

Opalustre Kit

ULT0554 

ENAMEL MICROABRASION



PREMA
PREMIER

PREMA microabrasion is a patented acid abrasion formulation used to correct superficial white, brown or multicoloured "demineralisation" enamel defects.

PREMA Kit: 10 cups, 2 mandrels, 10 mandrel tips, 1 finger ring.

PRE8009 



Opalescence PF
ULTRADENT



Whiten with Caries Prevention.

It's hard to beat science! Exciting findings of the Opalescence PF product line presented at the American Association of Dental Research (AADR) in 2006 confirm that whitening with Opalescence PF will provide ANTI-Caries benefits, Increase Enamel Microhardness, Minimise sensitivity and Improve Enamel Health.

Opalescence PF is highly recommended for whitening discoloured teeth prior to placement of composite, veneers, and/or crowns. Opalescence PF is effective in removing some or all of the internal tooth discolourations due to congenital, systemic, pharmacologic, traumatic, etc., factors and aging. Great success with fluorosis and even tetracycline staining.

The Opalescence PF sticky formula does not leach from the tray like other bleaching gels. Patients experience results quickly, which motivates them to pursue the treatment. Vital teeth whitened with Opalescence remain lighter for years. Opalescence PF is available in 4 concentrations, 4 formulations, 3 flavours and 2 kit configurations to meet all your whitening needs.

TOOTH WHITENING TAKE HOME



REV!
PREMIER

15 Minutes once a day

REV! whitening treatment is a revolutionary process that delivers optimal whitening in only 15 minutes, once a day.

Whitening Refresher Packs x 6
PRE0141 


Whitening Refresher Packs x 6
PRE0142 




Bravo
PREMIER

30 Minutes once a day

Shorter wear time only 7 hours is required to complete the 2 week take home treatment. No sensitivity patented water-based 9% hydrogen peroxide gel formulation hydrates teeth as they whiten to avoid sensitivity. Cosmetic patient pouch reflects well on your practice, and organises educational material allowing you to send home other practice building items.

Patient Pack 4x3ml Syringe
PRE0092 

Value Pack 24 x 3 ml Syringe, 6 Travel Pouches and 6 patient guides.
PRE0093 

Touch-Up Pack 1 x 3ml Syringe, dispensing tips.
PRE0091 

TOOTH WHITENING TAKE HOME

Doctor's Kits
ULTRADENT

An impressive presentation to patients and ideal for surgeries who purchase tray material in bulk or have a lab fabricate their trays

Doctors Kit contents:

- 8 x 1.2 ml prefilled syringes
- 1 x Tray case
- 1 oz Opalescence toothpaste
- 1.2 ml LC Block out
- 2 x 0.035" 5"x 5" Sof-tray sheets
- 1 x Shade Guide

- 10% Regular
ULT6509
- 10% Mint
ULT6510
- 10% Melon
ULT6511
- 15% Regular
ULT6512
- 15% Mint
ULT6513
- 15% Melon
ULT6514
- 20% Regular
ULT6515
- 20% Mint
ULT6516
- 20% Melon
ULT6517
- 35% Regular
ULT6201
- 35% Mint
ULT6204
- 35% Melon
ULT6207

Refill Kits
ULTRADENT

Suitable for surgeries that purchase materials in bulk and assemble their own kits.

- Contents 40 x 1.2ml prefilled syringe
- 10% Regular
ULT1306
 - 10% Mint
ULT1305
 - 10% Melon
ULT1307
 - 15% Regular
ULT1068

- 15% Mint
ULT1070
- 15% Melon
ULT1066
- 20% Regular
ULT1069
- 20% Mint
ULT1071
- 20% Melon
ULT1067
- 35% Regular
ULT6202
- 35% Mint
ULT6205
- 35% Melon
ULT6208



Ultra EZ® Desensitiser
ULTRADENT

The end to sensitivity!

UltraEZ® is ideal to use to quickly eliminate sensitivity from toothbrush abrasion, thermal and chemical changes. As well as tooth bleaching and root exposure. It is available in syringe to use with custom trays and disposable pre-loaded trays. UltraEZ® contains Potassium Nitrate and Fluoride to address both short and long term sensitivity.

- UltraEZ Syringe Kit
ULT1008
- UltraEZ Tray Delivered Upper
ULT7200
- UltraEZ Tray Delivered Lower
ULT7201



Tréswhite® Supreme
Ultradent

Pure 100% Opalescence PF gel in the convenient, Tréswhite delivery system.

New Opalescence Tréswhite Supreme is the first pre-loaded, convenient tray delivered over-the-counter take home whitening system in the market. Opalescence Tréswhite Supreme contains 10% hydrogen peroxide for faster and superior whitening results. It also contains potassium nitrate and fluoride. The convenience of pre-loaded trays of Opalescence Tréswhite Supreme means no lab and tray fabrication is required. This product is also ideal for patients where impression taking is a challenge. The sleek design of the trays provides ease of use and comfort to the patient. The external and internal bleaching trays can easily be positioned and adapted to fit any smile. Opalescence Tréswhite Supreme offers 3 delicious, improved flavours: mint, melon and peach. The optimum gel quantity of this product allows for rapid, easy cleanup following bleaching. The new high-tech packaging of Tréswhite provides superior stability and shelf life. Opalescence Tréswhite Supreme PF outstanding unique formula will give patients fresher and whiter teeth while rebuilding enamel and preventing caries. All Patient Kits Contain 10 upper and 10 lower pre-filled trays.

- Value Pack Mint - Pack of 6 Kits
ULT5706-Value Pack
- Value Pack Melon - Pack of 6 Kits
ULT5707-Value Pack
- Value Pack Peach - Pack of 6 Kits
ULT5708-Value Pack
- Patient Kit Mint
ULT5706
- Patient Kit Melon
ULT5707
- Patient Kit Peach
ULT5708

OPALESCENCE® MARKETING MATERIAL*

A full range of in surgery marketing material is available. Items shown are for illustration purposes. Contact your Gunz representative for a sample.



Posters

Made to fit frames, they make a great statement when displayed in surgeries and waiting rooms.



Waiting Room Brochures

Offer a quick and informative Q & A format on Tooth Whitening as well as before and after photographs.



Reception Stands & Mobiles

Designed to encourage questions from your patients.



Recall Cards

Send as post card appointment reminders to your patients.



Ceiling Tiles

Place these on the ceiling above dental chairs. Designed to excite patient interest in whitening procedures.

* Actual artwork may vary subject to availability.