



**GUNZDENTAL**

The people behind the products you trust

# NTI Heroes *Herald*

www.gunz.com.au

May/June 2007

## Extra! Extra!

# 25% DISCOUNT on NTI 'hero' products!\*



*Read all about it...*

\* Offer valid until 30 June, 2007. Discount only applies to NTI products featured in this brochure.

# Plaster

You can't go past these popular NTI instruments when bulk reduction or cutting through plaster is required.

CODE	PRODUCT	MAX RPM FOR INSTRUMENT	RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-365524450	<b>NTI Plaster cut diamond disc double sided 3mm</b>	10,000-15,000	\$92.84	<b>\$69.63</b>
NTI-365524300	<b>NTI Plaster cut diamond disc double sided 4.5mm</b>	10,000-15,000	\$61.05	<b>\$45.79</b>
NTI-HF251XCE060	<b>NTI Supercoarse Cross Cut size 060</b>		\$41.58	<b>\$31.19</b>
NTI-HF251GE060	<b>Bulk Cross Cutter size 060</b>		\$46.86	<b>\$35.15</b>

# Acrylic Grinders

CODE	PRODUCT	MAX RPM FOR INSTRUMENT	RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-HF108G060	<b>NTI Repair Cutter size 060</b> For roughening denture and acrylic teeth during repairs and remakes	15,000	\$44.11	<b>\$33.08</b>
NTI-AG430-075SC-HP	<b>NTI Acrylic Grinder AG430 HP 075SC</b> This instrument is perfect for smoothing and rounding the peripheries of custom trays, preventing soft tissue injuries to patients.	15,000	\$86.90	<b>\$65.18</b>
NTI-490544090	<b>NTI Diacryl grinder, 090, HP, Pkt 1</b> These acrylic trimmers help the technician to create adequate space in the denture for successful reline impressions.	10,000-15,000	\$50.93	<b>\$38.20</b>
NTI-P2100	<b>NTI Soft Wizard wheels x10</b> The only grinder you will need for all soft reline materials, custom bleaching trays and mouthguards.	15,000	\$43.12	<b>\$32.34</b>

# Acrylic Polishers

CODE	PRODUCT	MAX RPM FOR INSTRUMENT	RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-P0674	<b>NTI Acrylic Pre Polisher Coarse (Green)</b> Three-step silicone polishing system for prosthodontic, orthodontic and soft acrylics. Use as a replacement for sandpaper in difficult-to-reach areas.	15,000	\$32.56	<b>\$24.42</b>
NTI-P0664	<b>NTI Acrylic Pre Polisher Fine (Grey)</b>	15,000	\$32.56	<b>\$24.42</b>
NTI-P0654	<b>NTI Acrylic Pre Polisher High Gloss (Yellow)</b>	7,000	\$32.56	<b>\$24.42</b>
NTI-P0634	<b>NTI AcrylicMaster Coarse (Dark Blue)</b> Two-Step stone polishing system for polishing prosthetic acrylics. Coarse grit provides fast reduction and the Fine grit will help reduce polishing time.	10,000-15,000	\$32.56	<b>\$24.42</b>
NTI-P0644	<b>NTI AcrylicMaster Fine (Light Blue)</b>	10,000-15,000	\$32.56	<b>\$24.42</b>

# Tungsten Carbide Millennium Cutters

The smooth dimpled cutting surface of NTI Millennium Cutters enable ridge-free, smooth removal of metal chips. These instruments have an extreme metal hardness ensures higher efficiency, cooler cutting and longer working life.

CODE	PRODUCT	MAX RPM FOR INSTRUMENT	RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-MC079CE040	<b>NTI Millennium TC Cutter Standard Cross Cut Round Taper size 040</b>	40,000	\$49.61	<b>\$37.21</b>
NTI-MC251CE060	<b>NTI Millennium TC Cutter Standard Cross Cut Round Taper size 060</b>	40,000	\$49.61	<b>\$37.21</b>
NTI-MCL251CE060	<b>NTI Millennium TC cutter Standard Cross Cut Pointed Taper size 060</b>	40,000	\$49.61	<b>\$37.21</b>
NTI-MC079FE040	<b>NTI Millennium TC Cutter Fine Cross Cut Round Taper size 040</b>	40,000	\$58.63	<b>\$43.97</b>
NTI-MC251FE060	<b>NTI Millennium TC Cutter Fine Cross Cut Round Taper size 060</b>	40,000	\$79.53	<b>\$59.65</b>
NTI-MCL251FE060	<b>NTI Millennium TC Cutter Fine Cross Cut Pointed Taper size 060</b>	40,000	\$76.56	<b>\$57.42</b>
NTI-MC261FE023	<b>NTI Millennium TC Cutter Fine Cross Cut Pointed Taper size 023</b>	40,000	\$44.11	<b>\$33.08</b>

















## NEW Super Fine Tungsten Carbide Millennium Cutters

The unique cutting edge geometry of the new Super Fine cutters reduces vibration of the instrument by 95%, making it safer & easier when working in tight or delicate areas. Working with such low vibration levels helps to prevent wrist injuries and avoids the development of tendonitis and tennis elbow. Super Fine Millennium Cutters can be used for all dental alloys, porcelain, composite and denture acrylics.

CODE	PRODUCT	MAX RPM FOR INSTRUMENT	RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-MC073SFE014	<b>NTI Millennium TC Cutter Flame Super Fine size 14</b>	40,000	\$47.08	<b>\$35.31</b>
NTI-MC079SFE045	<b>NTI Millennium TC Cutter Round Taper Super Fine size 045</b>	30,000	\$79.97	<b>\$59.98</b>
NTI-MC129SFE023	<b>NTI Millennium TC Cutter Flat Taper Super Fine size 023</b>	40,000	\$49.61	<b>\$37.21</b>
NTI-MC138SFE008	<b>NTI Millennium TC Cutter Pointed Taper Super Fine size 008</b>	40,000	\$79.97	<b>\$59.98</b>
NTI-MC138SFE016	<b>NTI Millennium TC Cutter Pointed Taper Super Fine size 016</b>	40,000	\$79.97	<b>\$59.98</b>
NTI-MC138SFE023	<b>NTI Millennium TC Cutter Pointed Taper Super Fine size 023</b>	40,000	\$79.97	<b>\$59.98</b>
NTI-MCL077SFE023	<b>NTI Millennium TC Cutter Pear Super Fine size 023</b>	40,000	\$49.61	<b>\$37.21</b>
NTI-MC139SFE023	<b>NTI Millennium TC Cutter Bullet Super Fine size 023</b>	40,000	\$49.61	<b>\$37.21</b>
NTI-MC251SFE040	<b>NTI Millennium TC Cutter Flame Super Fine size 040</b>	30,000	\$79.97	<b>\$59.98</b>
NTI-MCL251SFE060	<b>NTI Millennium TC Cutter Flame Super Fine size 060</b>	20,000	\$79.97	<b>\$59.98</b>




# CeraGlaze Polishers

CeraGlaze is a 3-step polishing system for use with all types of porcelain. Unique binding material allows for optimal polishing by natural diamonds without the need for polishing paste. Minimal heat generation means CeraGlaze polishers can also be used on non-precious alloys and titanium. Highly recommended for CEREC and CAD/CAM restorations. For best results, use with light working pressure.

CODE	PRODUCT	MAX RPM FOR INSTRUMENT		RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-P310	CeraGlaze Polisher Green Pre Polish Slim Wheel	16,000		\$64.46	\$48.35
NTI-P342	CeraGlaze Polisher Green Pre Polish Disc	16,000		\$64.46	\$48.35
NTI-P344	CeraGlaze Polisher Green Pre Polish Point	16,000		\$64.46	\$48.35
NTI-P301	CeraGlaze Polisher Green Pre Polish Med Wheel	16,000		\$64.46	\$48.35
NTI-P3010	CeraGlaze Polisher Blue Refined Polish Slim Wheel	12,000		\$64.46	\$48.35
NTI-P3041	CeraGlaze Polisher Blue Refined Polish Bullet	12,000		\$64.46	\$48.35
NTI-P3042	CeraGlaze Polisher Blue Refined Polish Disc	12,000		\$64.46	\$48.35
NTI-P3043	CeraGlaze Polisher Blue Refined Polish Small Wheel	12,000		\$64.46	\$48.35
NTI-P3044	CeraGlaze Polisher Blue Refined Polish Point	12,000		\$64.46	\$48.35
NTI-P3001	CeraGlaze Polisher Blue Refined Polish Med Wheel	12,000		\$64.46	\$48.35
NTI-P30001	CeraGlaze Polisher Yellow High Shine Med Wheel	6,000		\$64.46	\$48.35
NTI-P30010	CeraGlaze Polisher Yellow High Shine Slim Wheel	6,000		\$64.46	\$48.35
NTI-P30041	CeraGlaze Polisher Yellow High Shine Bullet	6,000		\$64.46	\$48.35
NTI-P30042	CeraGlaze Polisher Yellow High Shine Disc	6,000		\$64.46	\$48.35
NTI-P30043	CeraGlaze Polisher Yellow High Shine Small Wheel	6,000		\$64.46	\$48.35
NTI-P30044	CeraGlaze Polisher Yellow High Shine Point	6,000		\$64.46	\$48.35




## All Ceramic SuperMax

The gentle trimming properties of the ceramic bonder reduces heat build-up making them perfect for trimming and smoothing delicate ceramic veneers.

CODE	PRODUCT	MAX RPM FOR INSTRUMENT		RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-G8001	NTI Super Max all ceramic Diam, Cylinder	10,000-15,000		\$38.72	\$29.04
NTI-G8002	NTI Super Max all ceramic Diam, Cone	10,000-15,000		\$38.72	\$29.04
NTI-G8003	NTI Super Max all ceramic Diam, Disc	10,000-15,000		\$57.86	\$43.40

## SuperMax Ceramic, Titanium, Chrome

High powered diamond grinders are cool cutting, self cleaning and self sharpening, making them the ideal instrument for titanium. Low dust generation and smooth cutting also is a great advantage when working with ceramics, chrome, composites and gold.

CODE	PRODUCT	MAX RPM FOR INSTRUMENT		RRP (incl. GST)	WITH 25% DISCOUNT (incl. GST)
NTI-G9001	NTI SUPER MAX Organ. bonded diamond Wheel	15,000-20,000		\$37.73	\$28.30
NTI-G9003	NTI SUPER MAX Organ. bonded diamond Cup	15,000-20,000		\$37.73	\$28.30
NTI-G9007	NTI SUPER MAX Organ. bonded diamond Disc	15,000-20,000		\$45.43	\$34.07

FOR ALL YOUR NTI NEEDS  
TALK TO US.



**GUNZDENTAL**  
The people behind the products you trust