



GUNZDENTAL

The people behind the products you trust

Products plus

ORDER FREEPHONE 1800 025 300 ORDER FREEFAX 1800 025 202 www.gunz.com.au

Issue 1 2005

Power-up!

Satelec® launches P5 Newtron generator

The latest in piezo ceramic technology

Already available on the Australian market place in an integrated format, Satelec's powerful "Newtron" technology will soon be available as a compact, portable benchtop unit.

Appropriately named, the **P5 Newtron**, offers all the benefits of the well known P5 Booster piezo generator, offering a wide range of applications from perio and endodontic applications, to general scaling, and, using the powerful booster range, heavy calculus removal and instrument retrieval.

The **Newtron technology** enables real-time tuning & adjustment (at 30,000 times per second) of the tips' single plane oscillation amplitude (Cruise Control™). This means less tip breakage and smooth and virtually painless vibrations for your patients, with unbeatable comfort during all treatments in prophylactics, endodontics and non-surgical periodontics.

Guided by Satelec's new **Colour Coding System** to select the required setting, the Newtron recognises each tip and calibrates automatically, guaranteeing optimal performance and reducing operator fatigue.

Along with the widest choice of ultrasonic instruments available in the marketplace today, Satelec also offers



the convenience of one Suprasson handpiece – interchangeable across the Satelec piezo generator range.

P5 Newtron is an easy-to-install, compact and reliable tabletop device. It has a smooth, modern aesthetic design, which is easy to clean and disinfect; and has been built in stringent compliance with the latest standards and safety regulations.

Satelec's Suprasson® P5 Newtron was previewed for the first time at IDS 2005 and will be available for sale from Gunz Dental in the second half of 2005.

The Newtron integrated piezo generator is currently available from Gunz Dental, order code: SATF55292

Call Gunz Dental on 1800 025 300 for further information.

sonicare
the sonic toothbrush

Now available in Australia




Philips has recently launched their Sonicare electric toothbrush in Australia. Sonicare is the brand personally used and recommended by more dental professionals in the USA and Japan.

Sonicare's patented technology creates a unique dynamic fluid cleaning action, making it exceptionally effective in hard-to-reach areas. Sonicare delivers a unique oral care experience with results you can see and feel. It is proven to remove more plaque biofilm beyond the bristles and to be 4 times gentler on dentin than oscillating/rotating toothbrushes.

Research results have confirmed many other interesting clinical facts about the Sonicare toothbrush.

For example:

- Patients brushed almost twice as long when using Sonicare as when using a manual toothbrush.



(cont. p2)

Inside:

P5 Booster now with Colour Coding	pg 2	Endo Eze AET. Responding to the Tooth	pg 6-7	Temporary repairs made easy with Luxatemp Fluorescence	pg 12
Sonicare launches in Australia		EndoRez is Revolutionising Root Canal Filling	pg 8	Get Optimal Occlusions with DMG O-Bite	pg 13
Introducing Satelec's new Endo & Scaling Tips	pg 3	Diamond Burs. Is Bigger Better?	pg 9	Streamlining your online ordering with ezyOrder	pg 14
Endo Twinn improving Obturation with Heat & Vibration	pg 4	Premier's unique Implant Cement for Long Term Temporaries	pg 10	Bordeaux is making more than Great Wine	pg 15
Tackling tough curves with Mani RT files	pg 5	Introducing Sword's latest implantology instruments	pg 11	News from Gunz	pg 16