

EveryDay EveryPatient

FEBRUARY - MARCH 2019



**SAVE
10%**

Medicom Chair Head Rest Covers

White tissue/poly. Easy to replace and remove repellent material.

MDC3014	25 x 25cm - 500 per carton	\$53.66
MDC3017	25 x 33cm - 500 per carton	\$61.13



**SAVE
10%**

Medicom SafeBasics™ Dry-Back® 2-Ply Bibs

Durable patient protection thanks to high absorption. Ideal for cleanings and minor procedures.

- Horizontal embossed pattern and a unique water-repellant edge for enhanced shielding from fluids.
- Strong, tear-resistant plastic backing allows for increased clothing protection.

MDC8282	Blue - Pack 500	\$37.92 per ctn
MDC8283	White - Pack 500	



**SAVE
15%**

Progress Nitrile White Gloves - Powder Free

Made from 100% Nitrile. Provide excellent protection with increased tensile strength and puncture resistance.

Carton - Contains 10 Boxes of Gloves	\$116.06 per ctn
--------------------------------------	-------------------------

Available Sizes: XS/S/M/L



**SAVE
15%**

Medicom Barrier Film Blue & Clear

Convenient dispensing with unique packaging that offers counter top dispenser with single use tear off sheets.

MDC5050	BLUE x1200 SHEETS	\$25.82
MDC5051	CLEAR x1200 SHEETS	

Clinell® Universal Wipes - The Professionals' Choice of Wipes!

Clinell Universal Wipes are one of the most effective antimicrobial products on the market. This revolutionary formula contains 6 different biocides. Each biocide has a different mechanism of action, reducing the risk of germs developing resistance.

Clinell Detergent Wipes provide an ideal choice for multi-surface general cleaning and damp dusting. Clinell Detergent Wipes are alcohol and disinfectant free.

AMCDW215	Clinell Neut Deterg Wipes Pack x215	\$14.39
AMCDT110	Clinell Neut Deterg Wipes 110 Tubx8	\$79.20
AMCW200	Clinell Sanitising Wipes 200pk	\$14.85
AMCWBU225	Clinell Sanitising Wipes Bucket 225	\$19.04
AMCWBU225R	Clinell Sanitising Wipes Bucket R/f	\$16.76



**SAVE
10%**





Medicom SafeSeal Quattro Sterilization Pouches

MDC88015	57x100mm - Pack 200
MDC88000	70x229mm - Pack 200
MDC88005	89x133mm - Pack 200
MDC88010	89x229mm - Pack 200
MDC88025	133x254mm - Pack 200
MDC88030	191x330mm - Pack 200
MDC88035	254x356mm - Pack 200
MDC88040	305x432mm - Pack 200

\$7.87
\$11.75
\$11.75
\$11.75
\$20.60
\$39.15
\$51.49
\$66.95



BUY 6 GET 1 FREE



3M Sterilisation Testing Consumable Refills

MMM1262P	3M Attest Biological Indicator Brown, Pkt of 25	\$223.50
MMM1300	3M Comply Disposable Bowie Dick Test Pack, Pkt 20	\$272.83

SAVE 5%

Steritech Bowie Dick Helix Test Set

MMP1012010100
 Steritech Bowie Dick Helix Test Indicators - 100
\$125.40

MMP1012010250
 Steritech Bowie Dick Helix Test Indicators - 250
\$245.58



SAVE 10%

Geringe Assured Bowie-Dick Test Card TableTop

Air Removal & Steam Penetration Test
 The Geringe Assured Bowie-Dick Test Card TT is for daily monitoring of small pre-vacuum type B steam sterilizers operating at 134°C. The card is specially packaged and designed to meet EN 867-5 requirements. This Geringe Assured product is tested and verified to have equivalent and approved performance at both EN 867-5 test cycles and Geringe Bowie-Dick cycles for type B pre-vacuum table top sterilizers.

AUSPL254 Bowie-Dick test Card TT - Pack 15 **\$48.84**



SAVE 10%

Geringe Titem Diagnostic Emulator

Class 6 Steam Strips
 Diagnostic emulator strip is designed to be used in steam autoclaves operating at 134 °C. They provide a precise integrated response to various defined combinations of temperature, time and the presence of steam.

AUSAUS508 Titem Emulating Steam Indicator Strips - Pack 250 **\$71.63**

Geringe Assured Helix Test – Kit

Heli-Checks – A hollow load challenge test for pre-vacuum steam sterilizers.

Hollow Load Testing:

- A good vacuum for steam penetration at proper time and temperature will change indicator colour from purple to green
- Reusable device for 250 uses
- Simple to load and remove indicator
- Package contains one device with resealable bag of 250 high-precision sterilant specific indicators
- 3-year shelf life
- Helix devices are not useful for testing of porous loads and therefore should not be used as a Bowie-Dick test

AUSLT104 Description **\$366.26**

SAVE 10%





**SAVE
10%**

Clinicare Neutral Detergent Towelettes

DL2610RZ	Refill - 220Wipes	\$12.95
DL2610R	Refills 200Wipes - Ctn 12	\$98.22
DL2610Z	Tub - 220Wipes	\$17.89
DL2610	Tubs 220Wipes - Ctn 12	\$112.58



**SAVE
10%**

Clinicare Hospital Grade Disinfectant Towelettes

DL2910RZ	Refill 220Wipes	\$19.37
DL2910R	Refills 220Wipes - Ctn 12	\$177.43
DL2910Z	Tub 220Wipes	\$23.57
DL2910	Tubs 220Wipes - Ctn 12	\$163.15



Handicare Handwash & Moisturiser

- Easy to dispense and set up in surgeries, clinics, hospitals and working environments either stand alone or using the simple bracket mount available
- Formulated for rapid, complete cleansing of the hands in preparation for clinical procedures and daily, repeated use
- pH 5.5 (natural pH of skin)
- Fragrance free – no residue fragrance or potential irritancy
- Colouring free – no residue colour, staining or potential irritancy
- Formula is completely biodegradable (environmentally safe)



**SAVE
15%**

DL1405	Moisturiser 500ml	\$11.71
DL1401	Handwash 500ml	\$10.22

V-Wipes Disinfectant Wipes

V-Wipes unique pH neutral formula and high quality apertured fabric allow for fast, effective cleaning and disinfection of all surfaces.

WHI210561
H.G Disinfectant wipes 100 wipes x4

\$62.33

WHI210579
V-Wipes Flat Packs 50 wipes x12

\$63.38



**SAVE
10%**

Whiteley Bactol Alcohol Gel

Bactol® Alcohol Gel provides antibacterial hand cleansing without the need for water.

WHI020066 with Macadamia Oil
500ml

\$14.90



**SAVE
10%**

Whiteley Hospital Grade Speedy Clean Wipes

Excellent ready-to-use hospital grade, anionic neutral detergent wipes which leave a streak free finish on surfaces. Ideal for cleaning and equipment in the surgery. The single-use application improves infection control by avoiding cross contamination.

WHI190328 1 x Canister of 100 **\$12.77**

WHI190333 1 x Flat Pack of 50 **\$4.33**



**SAVE
10%**



SAVE 15%

Mr Thirsty One-Step Isolation Device

Bring an inexpensive and efficient product like Mr. Thirsty One-Step into your dental practice and experience a hands free isolation device that retracts, isolates, and evacuates!

- IRC50Z987 Mr Thirsty One-Step Isolation Device Trial Kit
 - IRC50Z98710 Mr Thirsty One-Step Isolation Device Trial Kit 10mm
 - IRC50Z98716 Mr Thirsty One-Step Isolation Device Trial Kit 16mm
- \$52.97**



Progress G82+ Lubricant

High quality, synthetic oil suitable for use with ALL brands of dental handpieces.

PRO-G82-PLUS 1x 200g can **\$21.84**



BUY 3 GET 1 FREE

Dürr FD 366 Sensitive Wipes

Ready-to-use disinfectant and cleaning wipes for delicate surfaces. Bactericidal, tuberculocidal, yeasticidal, limited virucidal activity and non-enveloped viruses (noro)

- Fast disinfectant effect
- Extra-large wipes
- Perfume-free

Active ingredients: Alcohols and quaternary ammonium compounds
Wipe size: 17,5 x 25 cm

- DURFD366-TUB 100 Wipes Tub **\$108.55**
- DURFD366-REFILL 100 Wipes Refill - Ctn of 4 **\$326.20**

SAVE 10%



Dürr FD 366 Sensitive Rapid Disinfectant

Usable solution for the disinfection and cleaning of especially sensitive areas.

- Bactericidal, tuberculocidal, yeasticidal, limited virucidal
- Quick disinfectant effect
- For treatment units, ideal in combination with FD 360 artificial leather cleaner
- Tested material quality - recommended by leading unit manufacturers

Active substances: Alcohols and QAV

DURFD366 2.5 Litres **\$94.83**

Premier Brite Shield Instrument Cleaner

Suitable for use as a static soak or ultrasonic cleaner. Removes stains from instruments. 800gm.

PRE11100 Instrument Cleaner 800g **\$140.23**

SAVE 10%



SAVE 10%

Ultrasoak Dual Purpose Detergent

A dual-purpose neutral detergent concentrate used as an immersion solution to help remove materials from instruments.

PRO905 Neutral Detergent 5 Litre **\$25.64**

