



CREATE IT.



Phone: 02 9935 6677
Phone Orders: 1800 025 300
Email: dentalsales@gunz.com.au
Postal address: Locked Bag 5000 Alexandria NSW 2015

BUY 2 GET 1 FREE

March-April 2019

TI-MAX Z95L
3 YEAR WARRANTY

Ti·Max Z SERIES

Discover titanium

- Titanium Body
- Scratch Resistant DURAGRIP
- Cellular Glass Optics
- Clean Head System
- Push Button Chuck
- Microfilter
- Lighter Weight
- Slimmer Body
- Three Year Warranty



Ti-Max Z95L

OPTIC
ORDER CODE: C1034

- 1:5 Increasing
- Quattro spray
 - For FG burs (ø1.6)
 - Max speed: 200,000 min⁻¹
 - DLC coating
 - Anti-Heat System
 - Ceramic Bearings



Ti-Max Z25L

OPTIC
ORDER CODE: C1038

- 1:1 Direct drive
- Single spray
 - For CA burs (ø2.35)
 - Max speed: 40,000 min⁻¹



Ti-Max Z24L

OPTIC
ORDER CODE: C1126001

- 1:1 Direct drive
- Single spray
 - For CA burs (ø2.35)
 - Miniature Head
 - Max speed: 40,000 min⁻¹

NSK DISCOUNT OFFER

M600WLED
ORDER CODE: P1096
FOR W&H COUPLINGS
RRP \$1,472.63

BUY ONE
GET 40% OFF
\$880 each

40% OFF

BUY TWO
GET 50% OFF
\$749 each

50% OFF



Not to be used in conjunction with any other NSK promotions.

All quoted prices include GST.

NSK Oceania Pty Ltd www.nskoceania.com.au

Product Information | Technical Support
Freephone: 1300 44 33 21 (Australia) Email: info@nskoceania.com.au