

Experience the Difference!



PremierAir™

Hygiene & Diagnostic Instruments

Better Ergonomics, Less Fatigue.

new!

Single-ended Diagnostic
Explorers & Probes

PremierAir™ is here!

PremierAir™ is an innovative line of hygiene and diagnostic instruments that are specially designed to reduce hand fatigue. Ergonomic and light-weight, each instrument combines a rugged polymer handle with a comfortable non-slip grip. Less pressure is needed to maintain control and that means your hands and wrists will feel the difference!

Made to the highest standards by master craftspeople in the U.S.A., each instrument undergoes extensive testing to ensure a long working life. Every time.



PremierAir™

Scalers and Curettes

Hygiene Instruments

Precision tips made with 440A stainless steel and SmartSharp® Technology keep edges sharp and well defined.



Scalers & Curettes - 15g



PremierAir Columbia 4R/4L • 1004651



PremierAir Columbia 13/14 • 1004652



PremierAir Gracey 1/2 • 1004660



PremierAir Gracey 3/4 • 1004661



PremierAir Gracey 5/6 • 1004662




PremierAir Gracey 7/8 • 1004663



PremierAir Gracey 11/12 • 1004665



PremierAir Gracey 13/14 • 1004666




PremierAir Barnhart 5/6 • 1004679




PremierAir McCall 13S/14S • 1004671




PremierAir Nebraska 128 • 1004680



PremierAir McCall 17S/18S • 1004672




PremierAir U15/30 • 1004692




PremierAir H6/H7 • 1004689




PremierAir 137 • 1004695




PremierAir U15/33 • 1004693



PremierAir H204SD • 1004697



PremierAir 204S • 1004696



PremierAir NV2 • 1004682



PremierAir NV1 • 1004681




PremierAir Jacquette 34/35 • 1004686



PremierAir NV4 • 1004684




PremierAir H5/33 • 1004688




PremierAir H204S • 1004687

Premier's innovation in sharpening technology allows you to start sharp and work smart:


- ✓ Easily maintain sharpness – proprietary process ensures a well defined edge you can easily re-sharpen as needed
- ✓ Apply less pressure and use fewer strokes – achieve a glass-like tooth surface
- ✓ SmartSharp® available on Premier Scalers and Curettes – at no additional cost




Premier® Instrument Edge



Leading Competitor Edge



A sharp edge ensures a smooth, precise cut with less hand fatigue.



An irregular edge is less efficient and can lead to more hand fatigue.

PremierAir™

Diagnostic Instruments

new!
Single-end

Single-ended Explorers



new! PremierAir Explorer SE 23 • 1003551




new! PremierAir Explorer SE 23XT • 1003552



new! PremierAir Explorer SE 6 • 1003553

Single-ended Probes



new! PremierAir Probe SE Marquis 12 CV • 1003633



new! PremierAir Probe SE Williams 12 CV • 1003634



new! PremierAir Probe SE UNC 15 CV • 1003635



new! PremierAir Probe SE WHO CV • 1003636



new! PremierAir Probe SE UNC 12 CV • 1003637




SE Probes and Explorers - 9.5g




Double-ended Explorers

Explorers have excellent tactile sensitivity and sharp points for easy detection of tooth texture and surface irregularities.

DE Probex and Explorers - 12g




PremierAir Explorer DE 5 • 1003560



PremierAir Explorer DE CH3 • 1003563



PremierAir Explorer DE OD 11/12 • 1003562



PremierAir Explorer DE 6/23 • 1003565



PremierAir Explorer DE 2 • 1003564



PremierAir Explorer DE TU 17/23 • 1003566

Double-ended Probex


Probes' recessed dark markings are easy to read and resistant to wear.




PremierAir Probex Marquis 12 CV • 1003600
(3-6-9-12)



PremierAir Probex Williams CV • 1003601
(1-2-3-5-7-8-9-10)



PremierAir Probex UNC 12 CV • 1003604
(1-2-3-4-5-6-7-8-9-10-11-12)



PremierAir Probex UNC 15 CV • 1003602
(1-2-3-4-5-6-7-8-9-10-11-12-13-14-15)



PremierAir Probex WHO • 1003603
(3.5-5.5-8.5-11.5)

premier® Inspired solutions
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com