

## TECHNICAL DATA SHEET Saliva Ejector

### Intended Use

Saliva ejectors are used in dentistry in order to remove saliva, blood or other debris from the mouth during dental procedures. This disposable device is connected to an active device, normally a pneumatic system suction circuit.

Characteristic	Specification	
Material Composition	<b>Cannula and Filtering Part</b>	PVC (Polyvinyl Chloride), Latex free, Phthalate free
	<b>Wire Insert</b>	Zinc-coated iron wire
Dimensions	<b>Saliva Ejector Length</b>	150 mm
	<b>Cannula External Diameter</b>	6,55 mm
	<b>Tip Length</b>	14,8 mm
	<b>Tip External Diameter</b>	9 mm
Description	Flexible plastic cannula with metal wire insert and fixed filtering part (tip). Cannula has a smooth surface and filtering part has 6 openings.	

### Visual



### Additional Information

Storage Condition	Country of Origin
Store in a cool dry place	Italy

### Appendix 1 List of Item Codes and Attributes

## Appendix 1 (List of items)

*AMD Medicom Inc.*

Item Code	Description	Packaging Information	Color	Specific Attribute	Treatment
700	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	Clear with Blue Tip	Unscented	Non-Sterile
701	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	Clear with White Tip	Unscented	Non-Sterile
703	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	White with White Tip	Unscented	Non-Sterile
710	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	Dark Blue with Blue Tip	Unscented	Non-Sterile
711	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	Pink with Pink Tip	Unscented	Non-Sterile
712	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	Lavender with Lavender Tip	Unscented	Non-Sterile
713	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	N/A	Mint Scented	Non-Sterile
714	Medicom SafeBasics Saliva Ejectors	105 Ejectors/Bag, 10 Bags/Case, 1050 Ejectors/Case	N/A	Strawberry Scented	Non-Sterile

## Appendix 1 (List of items)

*AMD Medicom Inc.*

Item Code	Description	Packaging Information	Color	Specific Attribute	Treatment
700-EX	Medicom SafeBasics Saliva Ejectors	100 Ejectors/Bag, 10 Bags/Case, 1000 Ejectors/Case	Clear with Blue Tip	Unscented	Non-Sterile
701-EX	Medicom SafeBasics Saliva Ejectors	100 Ejectors/Bag, 10 Bags/Case, 1000 Ejectors/Case	Clear with White Tip	Unscented	Non-Sterile
703-EX	Medicom SafeBasics Saliva Ejectors	100 Ejectors/Bag, 10 Bags/Case, 1000 Ejectors/Case	White with White Tip	Unscented	Non-Sterile