Benda[®] Brush. The original, disposable applicator brush.

World leader in disposable brushes and preferred choice of dental offices since 1989.

Centrix revolutionized the disposable applicator market 16 years ago with its introduction of the Benda Brush, making it easier to apply whitening agents, etchants, bonding agents, sealants, desensitizers, adhesives and resins.

The brush's chemical resistant fibers hold up to the harshest dental materials and they allow for greater carrying capacity and better smoothing and spreading action.

Our brush is the only brush designed by a dentist for use in a dental office. At Centrix, we make dentistry easier.

- Molded head design bends to any angle with no bounce back
- Non-linting, non-absorbent, chemical resistant bristles
- Disposable, to eliminate clean-ups and cross contamination
- Precision manufactured in the USA for consistent quality
- Available in regular and mini head sizes
- · Variety of colors available for color coding procedures
- Single-ended or double-ended styles, with patented "snap apart" feature

Ordering is easy. Call 800-235-5862 or contact your preferred dealer. Visit centrixbendabrush.com for more information.







CE





ORDERING INFORMATION:

SINGLE END STYLE

BENDA BRUSH REGULAR

Single-ended bristle brush, regular point, black bristles REF 370101 1 x 144 Assorted colors

BENDA BRUSH BIG BOX

 REF 370501
 1 x 576 Red

 REF 370502
 1 x 576 Blue

 REF 370503
 1 x 576 Yellow

 REF 370504
 1 x 576 Green

BENDA BRUSH MINI Single-ended bristle brush, fine point,

white bristles REF 370111 1 x 144 Assorted colors

 BENDA BRUSH MINI BIG BOX

 REF 370601
 1 x 576 Red

 REF 370602
 1 x 576 Blue

 REF 370603
 1 x 576 Yellow

 REF 370604
 1 x 576 Green

Benda Brush, the world leader in disposable brushes since 1989.



DOUBLE END STYLE

BENDA TWIN™ Double-ended bristle brush, regular point, snap-apart handle REF 370105 200 (2 x 100) Pink REF 370106 200 (2 x 100) Blue REF 370107 400 (2 x 200) Refill, Pink REF 370108 400 (2 x 200) Refill, Blue

SUGGESTED ACCESSORIES

BONDLOGIC REF 310110 72 wells, 12 Black covers REF 310120 60 Black covers


Ordering is easy.

Call 1-800-235-5862 or contact your preferred dealer.

centrixdental.com







@ 2006 Centrix, Inc. All rights reserved. Benda^ $^{\otimes}$ is a registered trademark of Centrix, Inc. MK0020 REV C