

**WHEELCHAIR**  
*Platform*  
**WITH CORRECTA UNIT**



Available through:



P: 1800 025 300  
E: [dentalsales@gunz.com.au](mailto:dentalsales@gunz.com.au)

# Wheelchair Platform

The biggest challenge in the treatment of disabled patients is to place the patient in the unit chair. With the Gunz Dental Wheelchair lifter, the disabled patient can easily be brought to the treatment position with a single button without getting up from their chair.



With the Multiaxis Planetary headrest, the position that the patient is comfortable with can be easily adjusted.



Seat Safety switches prevent foot jam when lowering seat

# WHEELCHAIR *Platform* WITH CORRECTA UNIT

On dentist side, control panel, which is produced by hygienic material, provides activating cuspidor, cup filler, operating light and commanding also platform positioning.

Certified by :

DIN EN ISO 9001:2008  
DIN EN ISO 13485:2003  
Medical Devices Directives 93-42

CE 1984



With security switches prevents foot from getting stuck.



Multi Axis Planet headrest with free movement capability allows the most convenient position for patients.

Esthetic and ergonomic Correcta dentist tablet has 5 exits. Each exit can be designed as dentist requested instrument ( scaler, air/elec. Micromotor, fiber optic etc). Each exit's air and water regulations are separate.

Dental Lights optional photocell or, light emitting diode (LED) and movable head systems



## OPTIONAL

Peristaltic Pump for surgical operations with serum holder.

OPTIONAL  
Speed and direction information on the LCD display



Aspirator's filter can be removed easily for cleaning



With one button, it is possible to have water directly from unit distilled water tank or from city water



Removable and easy cleanable magnetic instruments pillow.



Silicon cover can be removed and completely autoclavable.







### **Multi axis movement**

Freedom in the working position  
Left right, in every angle  
possibility to fix the headrest.



### **Double Movement Headrest**

With double movement dentist can adjust easily  
most comfortable position for patient.  
For children adjust down,  
for wheelchair patient easy in reverse position.





**Soft upholstery option**

**What is your COLOR of life?**

Wide range selection of colors suitable for your working environment

