

New!

**Thoroughly
Cleans Surfaces**

**Promotes Strong, Reliable
Bonding**



UniClean™

Universal Restorative Cleaner

UniClean™ effectively removes salivary contamination from the intaglio surface after try-in – a critical step for achieving strong, long-lasting bonds. Its advanced formulation binds salivary proteins and creates a rinseable precipitate, ensuring a thoroughly clean surface.

- ✓ Compatible with all restorative materials: zirconia, glass ceramics, gold, metal, & composite
- ✓ Cleaning with UniClean™ promotes strong, reliable bonding performance – especially with zirconia.
- ✓ Fast 20-second application
- ✓ Ready to go from the bottle – no shaking required

- ✓ Optimized viscosity for accurate placement & easy handling

**Orange tint
ensures visibility!**

- ✓ For coverage & thorough rinsing



Available through your authorized dealer.

premier® Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • premierdentalco.com

©2025 Premier Dental | All Rights Reserved



UniClean™

Universal Restorative Cleaner


9007100 UniClean™

Contents: 5g bottle




*Your essential step for
after try-in cleaning!*

**Phosphate Contamination
on Zirconia vs. Cleaner***



**SBS to Saliva Contaminated Zirconia
After Cleaning With Different Cleaners***



ZR-Cem®

Self-Adhesive Universal Resin Cement

**Latest Technology =
Better Patient Experience**

- ✓ Innovative MDP formula enhances adhesion & bond to zirconia substrates
- ✓ Lasting color stability
- ✓ Virtually no post-operative sensitivity ensures patient comfort

ZR-Cem Standard Pack

One 8g Syringe, 10 Mix Tips,
10 Intraoral Tips

3001560 Shade White
 3001561 Shade A2
 3001562 Shade Translucent

ZR-Cem Refill Tips

3001418 35 mix tips & 35 intraoral tips



Premier® Universal Primer

Zirconia and Glass Ceramics, Metal and Composite Primer

Premier® Universal Primer, is a single-component adhesive primer that increases the retention between restorative substrates and resin cements.

- ✓ Ideal retention for zirconia, alumina, metal, glass ceramics + composite materials
- ✓ Dual coupling agents provides versatility and work with all substrates
- ✓ Easy to use – just apply, wait 60 seconds, and dry



3001568
5ml Bottle

premier® Inspired solutions
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

*Data on file.

Lit# 8610437 0825005 Rev1 NKPSM