



Dr's Light **PRIME**

Slim

&

Wide

Dual Wave

Next Level Tech  
& Revolutionary Design



## Slim & Wide

### Next LevelTech&Revolutionary Design




#### 5.2mm slim head height

It is easily accessible anywhere in the oral cavity.



#### Wide irradiating surface with a diameter of 12mm

Polymerization is possible even on a wide posterior surface at once.



#### Fast polymerization

Fast polymerization is possible with strong output.



#### Dual wavelength LED

Can polymerize for any kinds of material with dual wavelength (peak 405 & 460nm).



#### 360° twistable head

Can approach any where in oral cavity.



#### Autoon/off

Grip the device to turn it on, and release it to turn it off.



#### Use of interchangeable soft lenses

If the irradiation surface is damaged, it can be easily replaced.




#### Big capacity

High-quality, Big-capacity battery applied.



#### Uniform(Homogeneous) light and strong output

Uniform light is irradiated using 21 light emitting elements.



#### Color Touch LCD screen

It is easy to use by using the touch function.



#### Various curing mode(Ortho mode)

Six program modes that can be operated by touch are built-in.

