


Elara 28L: Largest capacity in a least space

Tray capacity



Elara11

	Tuttnauer Elara 11	Melag Vacuklav 43B+	Melag Vacuklav 23B+
Dimensions WxHxD	53 x 43 x 59.2 cm	46 x 50 x 68 cm	42.5 x 48.5 x 75 cm
Total volume	1.34 m ³	1.56 m ³	1.54 m ³
Chamber size Diameter; length	Ø 28cm ; length. 43 cm	Ø 25cm ; length. 45 cm	Ø 25cm ; length. 45 cm
Capacity	28.5 L	23 L	22 L
Trays	6	6	5 max.

Size comparison Elara 28L with Melag 23L and 22L:

With a **largest capacity**, Tuttnauer Elara 11 fits on a 60 cm bench with the **lowest volume space requirements**.

With **Elara 11**, you optimize your sterilization room and you get **extra capacity (28.5L)** compared to other autoclaves that offer between 17 and 23 litre chambers

